


i-CAB

i-CAB – RFID cabinet with real-time stock management & traceability



Do you want to manage the stock of your expensive consumables in real-time on item level with traceability up to the patient? Yes, you can!

Our **i-CAB** combined with our RFID TAGS glued onto the products and our **iOSS** communication software traces in real-time when a product is removed from the **i-CAB**, when products are supplied or when unused products are returned into the **i-CAB**.

In addition to this, the built-in **Access Control Module** only allows **i-CAB** access to authorised people. Our **i-BIN** on the other hand is doing the traceability of the products that are used for a patient during operation.

The **i-CAB** is functioning with all type of products, but is most suitable for active and real-time stock management of expensive products like implants, stents, etc....

The **i-CAB** is equipped with ISO 6040 **MODUL-iT** modules for optimal storage and can be interfaced with the central stock system through our **iOSS** communication software.

In case your products aren't RFID tagged by default, **matteo** delivers and installs the **i-CAB** with an RFID print module for RFID TAGS. These RFID TAGS can be glued onto the products.

Thanks to our i-CAB, active and real-time stock management becomes a piece of cake.

Advantages

- Real-time stock management on item level (incl. product info, lot and serial number, exp date)
- Real-time tracking of who supplies the **i-CAB** and who takes out/returns products
- All product movements are connected to the right person
- Product movements are shown on the **i-CAB**'s screen to ensure the registration was correct
- Real-time adjustment of the **i-CAB**'s stock level when products are taken, returned or supplied
- The pharmacist can follow the stock in the **i-CAB** in real time from the pharmacy
- The stock of consignment goods can be followed from distance by the external supplier
- Real-time traceability of used consumables per patient by throwing the packaging in our **i-BIN**
- Easy interface with the patient file software and the traceability software
- Save time thanks to real-time stock management and product traceability up to the patient.
- Easy procedures result in less mistakes and a better quality overall

Requirements per i-CAB

- Power plug AC 230 Volt
- Network connection RJ45 with fixed IP address

i-CAB – RFID cabinet with real-time stock management & traceability





mattteo

www.mattteo.com